



LifeGuard SmartTracker V2 (4G)

You want to remain independent, free to go about your normal daily activities and do the things you love. But your family worries when you're on your own. The solution? INS LifeGuard's newest Mobile Alarm System known as the SmartTracker 4G. With the simple press of a button, you're in contact with an experienced nurse. ***That's unique to INS LifeGuard!***

Simply press a button to send an alarm to a nurse in our 24/7 Emergency Response Centre. The nurse can speak with you via speakerphone and assess the situation. The built-in GPS locator ensures we'll be able to locate you in an emergency and send help — even if you are not able to respond.

The LifeGuard SmartTracker is perfect for seniors who like to be independent and active, as well as for WHS, lone workers or children. It is a stylish pendant that works using the same cellular phone technology as a mobile phone, but without the complicated features. It also functions as a fall detector, and has a docking station for easy charging.

Unlike other systems that depend on a friend or family member being available when you need help, our systems can send an alarm through to a nurse in our 24/7 Emergency Response Centre — *so help is always available.*

The LifeGuard SmartTracker provides assurance that you can contact a nurse at any time of the day or night — wherever you go. And your loved ones will have peace of mind knowing that help is at hand if you need it.



INS LifeGuard

*Caring for
You*

1800 636 226
www.theinsgroup.com.au

Features:

- Fall Detector*
- Set-up Geofences and Alerts that trigger if the device enters or exits the specified area
- Real-time location tracking - via 4 technologies: GPS, WiFi, Bluetooth and LBS
- A-GPS Support: Assisted GPS works using a network of satellites to quickly identify your location - faster than convention GPS, so help can arrive faster
- Bluetooth 5.0 - connect to your charging base for indoor locating

* Not all falls can be detected, therefore you should always push the HELP button after a fall, if you're able to.

Specifications:

4G LTE	GPS / GLNSS
<ul style="list-style-type: none">• Supported Bands:<ul style="list-style-type: none">» AU/SA: LTE B1/B3/B5/B7/B28; WCDMA B1/B5» EU: B1/B3/B7/B8/B20; WCDMA B1/B8; GSM 900/1800» North America variants: B2/B4/B5/B12/B17/B13; WCDMA B2/B4/B5• Antenna: Built-in laser curved antenna	<ul style="list-style-type: none">• GPS chipset: MTK6625 (AGPS support)• Support: GPS and Glonass• Antenna: Built-in ceramic antenna

General Specifications	Hardware Specifications
<ul style="list-style-type: none">• Dimensions: 62mm * 47mm * 18mm• Weight: 54g• Backup battery: Rechargeable, 3.7V, 1000mAh• Charging voltage: 5V DC• Operating temperature: -20°C to +80°C for working, -30°C to +70°C for storage• Battery life: Up to 120 hours under normal environmental conditions• Water resistant: Splashproof (not to be submerged or worn in the shower)• Warranty: 90 days	<ul style="list-style-type: none">• MCU: MTK6739, Nordic NRF52840• Sensors: Motion and Vibration Sensors• Connectors: 4 Pin-Magnet for Charging• SIM Card Slot: Nano SIM Card• Flash Memory: Built-in 1 MB memory• Built-in microphone and speaker• WiFi: 802.11 b/g/n, 2.4G• BLE: BT5.0 LE

Specifications subject to change at discretion of INS LifeGuard.

Please Note:

- Whilst the SmartTracker may work indoors, it's designed for outdoor use as it can be difficult to get a GPS signal indoors.
- Whilst every care has been taken to ensure our products and services meet high quality standards, the Lifeguard SmartTracker relies on mobile network connectivity at the time an alarm is being sent. As such, we cannot provide an absolute guarantee that it will function as expected at all times
- INS does not monitor battery, network or hardware issues.

Warning for those with pacemakers:

Radiofrequency energy from mobile devices and other technologies, including our transmitters and the LifeGuard SmartTracker, may cause interference with pacemakers. You should consult with your doctor to be sure, however most who wear a pacemaker can still safely use these devices by taking these simple precautions:

- **PROPER STORAGE** – avoid placing or wearing them next to the pacemaker. For instance, don't carry your mobile phone in your shirt or jacket pocket if it rests near the pacemaker. Wear your transmitter on the wrist rather than as a pendant.
- **PROPER DISTANCE** – hold your mobile phone to the ear opposite the side of the body where the pacemaker is implanted. Ensure your transmitters and/or the LifeGuard SmartTracker remain at least 15 centimetres away from your pacemaker.